DRUG INCORPORATION OF ACRYLIC RESINS MICROBIOLOGICAL AND RELEASE STUDIES

FARMÁCIA



COSTA J.¹, ALEXANDRE F.¹, MARCELINO N.¹, RIBEIRO I.², BETTENCOURT A.², PORTUGAL J.¹, NEVES CB.¹

¹ Biomedical and Oral Sciences Research Unit (UICOB), Faculty of Dental Medicine, University of Lisbon, Lisbon, Portugal ²Research Institute for Medicines (iMED.ULisboa), Faculty of Pharmacy, University of Lisbon, Lisbon, Portugal

PURPOSE

Denture stomatitis is a pathological condition of the denture bearing mucosa with multifactorial etiology, being the infection by Candida species especially Candida albicans considered the main etiologic factor¹⁻³. Inhibition of the formation of *C. albicans* biofilms on prostheses may be very important in preventing the development of denture stomatitis; therefore, chlorhexidine (CHX) has been incorporated in acrylic resins^{4,5}. Currently, only mechanical properties of acrylic reline resins incorporated with CHX were studied^{2,6,7}.

OBJECTIVES

Evaluate the drug release and antifungal activity against Candida albicans of acrylic reline resins loaded with chlorhexidine.

MATERIALS AND METHODS

ACRYLIC RESINS








PREPARATION OF **SPECIMENS**



CHX Diacetate

0%, 1%, 2.5%, 5%, 7.5% e 10% (w/w)




DRUG RELEASE





 Releasing solutions Artificial saliva











UV-spectroscopy (255nm)

ANTIFUNGAL ACTIVITY — Kirby-Bauer test

1) Inoculum Preparation: Direct Colony Suspension Method












2) Inoculation of Test Plates







3) Application of disks to inoculated agar plates




EXPERIMENTAL GROUPS C – Resin with 0% CHX 1 – Resin with 1% CHX 2,5 - Resin with 2.5% CHX 5 – Resin with 5% CHX 7,5 – Resin with 7.5% CHX 10 – Resin with 10% CHX

POSITIVE CONTROLS CHX – Paper disk with 10µg CHX F – Paper disk with 20µg Amphotericin B

Incubation at 35 ±2 °C for 48 hours


4) Reading plates Measurement of the diameters of inhibition zones (mm)





RESULTS


DRUG RELEASE



For all CHX % and for the majority of time intervals, Ufi Gel Hard released the highest amount of CHX, followed by Kooliner and **Probase Cold.**

ANTIFUNGAL ACTIVITY

No group with 1% CHX presented inhibition zone, and only the **Kooliner** group presented inhibition zone in specimens with 2.5% CHX (8.95±0.97 mm).



Probase Cold presented lower values than Ufi Gel Hard (p=0.017).

There was a tendency to increase the diameter of inhibition zone with the increase of the percentage of CHX incorporated in the material (p=0.198).

CONCLUSION

- Different acrylic reline resins compositions and different CHX loading percentages affects the drug release.
- The incorporation of CHX into acrylic resins seems to have an influence on the microbiological activity against this strain of C. albicans.
- U appears to be the resin which releases more CHX and which, at 5% CHX concentration, already exhibits antifungal activity.
- CHX delivery systems based on acrylic reline resins are a potential approach in the treatment of denture stomatitis.

vitro analysis of antifungal activity, drug release and hardness. Dental materials: official publication of the Academy of Dental Materials: official publication of the Academy of Dental Materials: official publication of the Academy of Dental Materials of incorporating antifungals on the water sorption and solubility of interim resilient liners

for denture base relining. The Journal of prosthetic dentistry. 2016;115(5):611-6. 7 - Urban VM, Lima TF, Bueno MG, Giannini M, Arioli Filho JN, de Almeida AL, et al. Effect of the addition of antimicrobial agents on Shore A hardness and roughness of soft lining materials. Journal of prosthodontics: official journal of the

American College of Prosthodontists. 2015;24(3):207-14. 8 - NCCLS. National Committee for Clinical Laboratory Standards. Publication of M44-A2: Method for Antifungal Disk Diffusion Susceptibility Testing of Yeasts: approved guideline. Second edition. NCCLS 2002;24(15):1-23.